

Product datasheet for **HP207352**

AP1S2 Human qPCR Primer Pair (NM_003916)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8905
Forward Sequence:	CAAGAGAACTTGTCAGACCGTTT
Reverse Sequence:	GTCCTGATCCTCAATAGCACAGC
ACCN:	NM_003916 , NM_003916.1 , NM_003916.2 , NM_003916.3 , NM_003916.4 , BC071867 , BC071867.1 , BC001117 , BI462192 , BT006738 , BX537780 , NM_003916.5
UniProt ID:	P56377
Synonyms:	DC22; MRX59; MRXS5; MRXS21; MRXSF; PGS; SIGMA1B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)