

Product datasheet for **HP207376**

CDKL2 Human qPCR Primer Pair (NM_003948)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8999
Forward Sequence:	GATTTGCGCGAACATTGGCAGC
Reverse Sequence:	TCAGTTACCAGACAACCAATGGC
ACCN:	<u>NM_003948</u> , <u>NM_003948.1</u> , <u>NM_003948.2</u> , <u>NM_003948.3</u> , <u>NM_003948.4</u> , <u>BC093646</u> , <u>BC093646.1</u> , <u>BC017417</u> , <u>BC093981</u> , <u>BX095447</u> , <u>NM_003948.5</u>
UniProt ID:	<u>Q92772</u>
Synonyms:	KKIAMRE; P56
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)