

Product datasheet for **HP207389**

CCRL2 Human qPCR Primer Pair (NM_003965)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9034
Forward Sequence:	TGCCGCTGTTCCATCTGCGTA
Reverse Sequence:	ACACTTCGGTGGAATGGTCAGG
ACCN:	<u>NM_003965</u> , <u>NM_003965.1</u> , <u>NM_003965.2</u> , <u>NM_003965.3</u> , <u>NM_003965.4</u> , <u>BC071682</u> , <u>BC096075</u> , <u>BC096076</u> , <u>NM_003965.5</u>
UniProt ID:	<u>O00421</u>
Synonyms:	ACKR5; CKRX; CRAM; CRAM-A; CRAM-B; HCR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)