

Product datasheet for **HP207409**

NFKB1 Human qPCR Primer Pair (NM_003998)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4790
Forward Sequence:	GCAGCACTACTTCTTGACCACC
Reverse Sequence:	TCTGCTCCTGAGCATTGACGTC
ACCN:	<u>NM_003998</u> , <u>NM_003998.1</u> , <u>NM_003998.2</u> , <u>NM_003998.3</u> , <u>BC051765</u> , <u>BC033210</u> , <u>BQ023914</u>
UniProt ID:	<u>P19838</u>
Synonyms:	CVID12; EBP-1; KBF1; NF-kappa-B1; NF-kappaB; NF-kappabeta; NF-kB; NF-kB1; NFkappaB; NFKB-p50; NFKB-p105
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)