

Product datasheet for **HP207410**

OSMR Human qPCR Primer Pair (NM_003999)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9180
Forward Sequence:	CAGGTGTTCTACCAAATCTGCG
Reverse Sequence:	AATCCACCCTCTGTGCCTGCAA
ACCN:	NM_003999 , NM_003999.1 , NM_003999.2 , BC125209 , BC010943 , BC063468 , BC125210 , BM893608 , BP352496 , BQ772952 , BU619548 , NM_003999.3
UniProt ID:	Q99650
Synonyms:	IL-31R-beta; IL-31RB; OSMRB; OSMRbeta; PLCA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)