

## Product datasheet for **HP207412**

### ACOX1 Human qPCR Primer Pair (NM\_004035)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51
Forward Sequence:	GGCGCATACATGAAGGAGACCT
Reverse Sequence:	AGGTGAAAGCCTTCAGTCCAGC
ACCN:	<a href="#">NM_004035</a> , <a href="#">NM_004035.1</a> , <a href="#">NM_004035.2</a> , <a href="#">NM_004035.3</a> , <a href="#">NM_004035.4</a> , <a href="#">NM_004035.5</a> , <a href="#">NM_004035.6</a> , <a href="#">BC008767</a> , <a href="#">BC008767.2</a> , <a href="#">BC010425</a> , <a href="#">BQ943492</a> , <a href="#">BX537380</a> , <a href="#">NM_004035.7</a>
UniProt ID:	<a href="#">Q15067</a>
Synonyms:	ACOX; MITCH; PALMCOX; SCOX
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)