

## Product datasheet for **HP207420**

### ATOX1 Human qPCR Primer Pair (NM\_004045)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	475
Forward Sequence:	TCTGAGCACAGCATGGACACTC
Reverse Sequence:	TCTGGAAGCCAGCGGGAGGAT
ACCN:	<a href="#">NM_004045</a> , <a href="#">NM_004045.1</a> , <a href="#">NM_004045.2</a> , <a href="#">NM_004045.3</a> , <a href="#">BC112248</a> , <a href="#">BC112250</a> , <a href="#">BI599233</a> , <a href="#">BT009786</a> , <a href="#">BU172675</a> , <a href="#">NM_004045.4</a>
UniProt ID:	<a href="#">O00244</a>
Synonyms:	ATX1; HAH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)