

## Product datasheet for **HP207423**

### BDH1 Human qPCR Primer Pair (NM\_004051)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	622
Forward Sequence:	GAAAGTGGTGGAGATTGTCCGC
Reverse Sequence:	TGTAGGTCTCCAGGCTGGTGAA
ACCN:	<a href="#">NM_004051</a> , <a href="#">NM_004051.1</a> , <a href="#">NM_004051.3</a> , <a href="#">NM_004051.4</a> , <a href="#">BC005844</a> , <a href="#">BC011964</a> , <a href="#">BC019317</a> , <a href="#">BC021183</a> , <a href="#">BE782251</a> , <a href="#">BI601294</a> , <a href="#">BQ897235</a> , <a href="#">NM_004051.5</a>
UniProt ID:	<a href="#">Q02338</a>
Synonyms:	BDH; SDR9C1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)