

## Product datasheet for **HP207432**

### CDH12 Human qPCR Primer Pair (NM\_004061)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1010
Forward Sequence:	CTGTCACCTGCTCAAGACCTGGA
Reverse Sequence:	GTGGCGATGGTTCCTTCATTTCC
ACCN:	<u><a href="#">NM_004061</a></u> , <u><a href="#">NM_004061.1</a></u> , <u><a href="#">NM_004061.2</a></u> , <u><a href="#">NM_004061.3</a></u> , <u><a href="#">NM_004061.4</a></u> , <u><a href="#">BC047608</a></u> , <u><a href="#">BC047608.1</a></u> , <u><a href="#">NM_004061.5</a></u>
UniProt ID:	<u><a href="#">P55289</a></u>
Synonyms:	CDHB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)