

Product datasheet for **HP207444**

Citrate synthetase (CS) Human qPCR Primer Pair (NM_004077)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1431
Forward Sequence:	CACAGGGTATCAGCCGAACCAA
Reverse Sequence:	CCAATACCGCTGCCTTCTCTGT
ACCN:	<u>NM_004077</u> , <u>NM_004077.1</u> , <u>NM_004077.2</u> , <u>BC072016</u> , <u>BC072016.1</u> , <u>BC000105</u> , <u>BC010106</u> , <u>BI258893</u> , <u>BI553413</u> , <u>BT007414</u> , <u>BX640838</u> , <u>NM_004077.3</u>
UniProt ID:	<u>O75390</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)