

Product datasheet for **HP207452**

TIMM8A Human qPCR Primer Pair (NM_004085)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1678
Forward Sequence:	GGTAGAGACTCAAAGCAGCGC
Reverse Sequence:	GAAGCGCTCAACGCAGTTCACA
ACCN:	NM_004085 , NM_004085.1 , NM_004085.2 , NM_004085.3 , BC015093 , BC005236 , BC006994 , BC070284 , BM467820 , BQ013177 , NM_004085.4
UniProt ID:	O60220
Synonyms:	DDP; DDP1; DFN1; MTS; TIM8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)