

Product datasheet for **HP207462**

STOM Human qPCR Primer Pair (NM_004099)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2040
Forward Sequence:	TCTGCCATGCACTGACAGCTTC
Reverse Sequence:	GAACGCGTAATAGACCACACC
ACCN:	NM_004099 , NM_004099.1 , NM_004099.2 , NM_004099.3 , NM_004099.4 , NM_004099.5 , BC010703 , BC010703.1 , BG701899 , BI553227 , BI763647 , BM451470 , BM925356 , BP208523
UniProt ID:	P27105
Synonyms:	BND7; EPB7; EPB72
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)