

Product datasheet for **HP207465**

FABP3 Human qPCR Primer Pair (NM_004102)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2170
Forward Sequence:	GTGGAGTTCGATGAGACAACAGC
Reverse Sequence:	TGGTCTCTGCCCGTCCCATTT
ACCN:	NM_004102 , NM_004102.1 , NM_004102.2 , NM_004102.3 , NM_004102.4 , BC007021 , BC007021.1 , BF576953 , BI752278 , BP194697 , BT006727 , NM_004102.5
UniProt ID:	P05413
Synonyms:	FABP11; H-FABP; M-FABP; MDGI; O-FABP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)