

Product datasheet for **HP207478**

FKBP1B Human qPCR Primer Pair (NM_004116)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2281
Forward Sequence:	GGCAAACAGGAAGTCATCAAAGG
Reverse Sequence:	TGGCTCCATATGCCACATCAGG
ACCN:	<u>NM_004116</u> , <u>NM_004116.1</u> , <u>NM_004116.2</u> , <u>NM_004116.3</u> , <u>NM_004116.4</u> , <u>BC002614</u> , <u>BC002614.1</u> , <u>BC050998</u> , <u>BF216711</u> , <u>BF969580</u> , <u>BI546013</u> , <u>BM802489</u> , <u>BU147063</u> , <u>BU508033</u> , <u>NM_004116.5</u>
UniProt ID:	<u>P68106</u>
Synonyms:	FKBP1L; FKBP12.6; OTK4; PKBP1L; PPlase
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)