

Product datasheet for **HP207489**

TFIIF (GTF2F2) Human qPCR Primer Pair (NM_004128)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2963
Forward Sequence:	GAAACTGCGGATTGCCAAGACTC
Reverse Sequence:	ATGTCTGTCTCCAACACTTTGC
ACCN:	<u>NM_004128</u> , <u>NM_004128.1</u> , <u>NM_004128.2</u> , <u>BC001771</u> , <u>BC001771.1</u> , <u>BT019525</u> , <u>NM_004128.3</u>
UniProt ID:	<u>P13984</u>
Synonyms:	BTF4; RAP30; TF2F2; TFIIF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)