

Product datasheet for **HP207492**

Granzyme B (GZMB) Human qPCR Primer Pair (NM_004131)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3002
Forward Sequence:	CGACAGTACCATTGAGTTGTGCG
Reverse Sequence:	TTCGTCCATAGGAGACAATGCC
ACCN:	NM_004131 , NM_004131.2 , NM_004131.3 , NM_004131.4 , NM_004131.5 , BC030195 , BC030195.1 , BX283601 , NM_004131.6
UniProt ID:	P10144
Synonyms:	C11; CCPI; CGL-1; CGL1; CSP-B; CSPB; CTLA1; CTSL1; HLP; SECT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)