

Product datasheet for **HP207501**

LLGL1 Human qPCR Primer Pair (NM_004140)

Product data:

| | |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 3996 |
| Forward Sequence: | CCTCTTCCAGACAGACTGTGAG |
| Reverse Sequence: | AGAGCAACCTTCTGCACGCCAA |
| ACCN: | <u>NM_004140</u> , <u>NM_004140.1</u> , <u>NM_004140.2</u> , <u>NM_004140.3</u> , <u>BC028037</u> , <u>BC028214</u> , <u>BC051466</u> , <u>BC151838</u> , <u>BC167831</u> , <u>BQ923585</u> , <u>BQ937224</u> |
| UniProt ID: | <u>Q15334</u> |
| Synonyms: | DLG4; HUGL; HUGL-1; HUGL1; Lgl1; LLGL; Mgl1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)