

Product datasheet for **HP207516**

RAB1A Human qPCR Primer Pair (NM_004161)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5861
Forward Sequence:	GGGAACAAATGTGATCTGACCAC
Reverse Sequence:	GAAAGACTGTTCTACATTGTTGC
ACCN:	<u>BC000905</u> , <u>NM_004161</u> , <u>NM_004161.1</u> , <u>NM_004161.2</u> , <u>NM_004161.3</u> , <u>NM_004161.4</u> , <u>BC000905.2</u> , <u>BG704216</u> , <u>BT019528</u> , <u>BT019529</u> , <u>BX571747</u> , <u>NM_004161.5</u>
UniProt ID:	<u>P62820</u>
Synonyms:	RAB1; YPT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)