

## Product datasheet for **HP207534**

### TANK Human qPCR Primer Pair (NM\_004180)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10010
Forward Sequence:	CCACTTCTGGACCCATCTGATG
Reverse Sequence:	GCAGTTCTGAGTCTGTGCCACT
ACCN:	<u>NM_004180</u> , <u>NM_004180.1</u> , <u>NM_004180.2</u> , <u>BC067779</u> , <u>BC067779.1</u> , <u>BC003388</u> , <u>BT009855</u> , <u>NM_004180.3</u>
UniProt ID:	<u>Q92844</u>
Synonyms:	I-TRAF; ITRAF; TRAF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)