

## Product datasheet for **HP207535**

### PGP9.5 (UCHL1) Human qPCR Primer Pair (NM\_004181)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7345
Forward Sequence:	CAGTTCAGAGGACACCCTGCTG
Reverse Sequence:	CCACAGAGCATTAGGCTGCCTT
ACCN:	<a href="#">NM_004181</a> , <a href="#">NM_004181.1</a> , <a href="#">NM_004181.2</a> , <a href="#">NM_004181.3</a> , <a href="#">NM_004181.4</a> , <a href="#">BC000332</a> , <a href="#">BC000332.2</a> , <a href="#">BC005117</a> , <a href="#">BC006305</a> , <a href="#">BC018838</a> , <a href="#">BP202166</a> , <a href="#">NM_004181.5</a>
UniProt ID:	<a href="#">P09936</a>
Synonyms:	HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)