

Product datasheet for **HP207541**

KDM5C Human qPCR Primer Pair (NM_004187)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8242
Forward Sequence:	ACTGCTGACCATTGCTGAACGC
Reverse Sequence:	CCTCCTTGAGAGCCTGGATGTT
ACCN:	<u>NM_004187</u> , <u>NM_004187.1</u> , <u>NM_004187.2</u> , <u>NM_004187.3</u> , <u>BC054499</u> , <u>BC047751</u> , <u>BI711169</u> , <u>NM_004187.5</u>
UniProt ID:	<u>P41229</u>
Synonyms:	DXS1272E; JARID1C; MRX13; MRXJ; MRXSCJ; MRXSJ; SMCX; XE169
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)