

Product datasheet for **HP207545**

ASMTL Human qPCR Primer Pair (NM_004192)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8623
Forward Sequence:	CTGCGGATATTGCCAGCAAAGTG
Reverse Sequence:	GAGTATTCGCCATCCGATGCCA
ACCN:	NM_004192 , NM_004192.1 , NM_004192.2 , NM_004192.3 , BC002508 , BC002508.2 , BC010089 , BM021926 , BM129217 , BU688748 , NM_004192.4
UniProt ID:	O95671
Synonyms:	ASMTLX; ASMTLY; ASTML
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)