

Product datasheet for **HP207581**

KLF4 Human qPCR Primer Pair (NM_004235)

Product data:

| | |
|--------------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 9314 |
| Forward Sequence: | CATCTCAAGGCACACCTGCGAA |
| Reverse Sequence: | TCGGTCGCATTTTGGCACTGG |
| ACCN: | <u>BC030811</u> , <u>NM_004235</u> , <u>NM_004235.1</u> , <u>NM_004235.2</u> , <u>NM_004235.3</u> , <u>NM_004235.4</u> , <u>NM_004235.5</u> , <u>BC030811.1</u> , <u>BC029923</u> , <u>BU739190</u> , <u>NM_004235.6</u> |
| UniProt ID: | <u>O43474</u> |
| Synonyms: | EZF; GKLF |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)