

## Product datasheet for **HP207586**

### Cip4 (TRIP10) Human qPCR Primer Pair (NM\_004240)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9322
Forward Sequence:	CGAATATGCGGCTCAACTGCAG
Reverse Sequence:	CCTGCGTTCATCCATGTCTTGG
ACCN:	<a href="#">NM_004240</a> , <a href="#">NM_004240.1</a> , <a href="#">NM_004240.2</a> , <a href="#">NM_004240.3</a> , <a href="#">BC013002</a> , <a href="#">BC013002.2</a> , <a href="#">BG494851</a> , <a href="#">BQ931996</a> , <a href="#">BT006698</a> , <a href="#">BT020167</a> , <a href="#">BT020171</a> , <a href="#">NM_004240.4</a>
UniProt ID:	<a href="#">Q15642</a>
Synonyms:	CIP4; HSTP; STOT; STP; TRIP-10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)