

## Product datasheet for **HP207588**

### **HMGN3 Human qPCR Primer Pair (NM\_004242)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	9324
<b>Forward Sequence:</b>	CACAAGACGGTCTGCCAGATTG
<b>Reverse Sequence:</b>	CCTCCTTCTCCCTTTAGCACC
<b>ACCN:</b>	<a href="#">NM_004242</a> , <a href="#">NM_004242.1</a> , <a href="#">NM_004242.2</a> , <a href="#">NM_004242.3</a> , <a href="#">BC009529</a> , <a href="#">BC009529.2</a> , <a href="#">BF790582</a> , <a href="#">BI551210</a> , <a href="#">BK000003</a> , <a href="#">BK000006</a> , <a href="#">BM927225</a> , <a href="#">BU619570</a> , <a href="#">BX648085</a> , <a href="#">NM_004242.4</a>
<b>UniProt ID:</b>	<a href="#">Q15651</a>
<b>Synonyms:</b>	PNAS-24; PNAS-25; TRIP7
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)