

## Product datasheet for **HP207602**

### SEP15 Human qPCR Primer Pair (NM\_004261)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9403
Forward Sequence:	TGTTGGCGACTGTGCTTCAAGC
Reverse Sequence:	GAAGCAGGTTGAACTGTCCGAG
ACCN:	<u><a href="#">NM_004261</a></u> , <u><a href="#">NM_004261.1</a></u> , <u><a href="#">NM_004261.2</a></u> , <u><a href="#">NM_004261.3</a></u> , <u><a href="#">NM_004261.4</a></u> , <u><a href="#">BC021697</a></u> , <u><a href="#">BC021697.1</a></u> , <u><a href="#">BC016359</a></u> , <u><a href="#">BC016359.1</a></u> , <u><a href="#">BC005294</a></u> , <u><a href="#">BQ429619</a></u> , <u><a href="#">BU941543</a></u> , <u><a href="#">NM_004261.5</a></u>
UniProt ID:	<u><a href="#">O60613</a></u>
Synonyms:	SEP15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)