

Product datasheet for **HP207612**

CHST3 Human qPCR Primer Pair (NM_004273)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9469
Forward Sequence:	GAAGCAAATACGCCCTTTTCTTGG
Reverse Sequence:	TCGGTGCTGTTGGCATCTGCTA
ACCN:	NM_004273 , NM_004273.1 , NM_004273.2 , NM_004273.3 , NM_004273.4 , BC093690 , BC093690.1 , BC033183 , BC104856 , BG325090 , NM_004273.5
UniProt ID:	Q7LGC8
Synonyms:	C6ST; C6ST1; HSD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)