

Product datasheet for **HP207613**

AKAP6 Human qPCR Primer Pair (NM_004274)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9472
Forward Sequence:	GCTCTGGATACCATCAAAGCCG
Reverse Sequence:	GAGGGACATCTGCTCAACAGCA
ACCN:	<u>NM_004274</u> , <u>NM_004274.1</u> , <u>NM_004274.2</u> , <u>NM_004274.3</u> , <u>NM_004274.4</u> , <u>BC137232</u> , <u>BC137233</u> , <u>BC150185</u> , <u>BC150288</u> , <u>BC154413</u> , <u>BP230135</u> , <u>NM_004274.5</u>
UniProt ID:	<u>Q13023</u>
Synonyms:	ADAP6; ADAP100; AKAP100; mAKAP; PRKA6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)