

Product datasheet for **HP207636**

ABCB7 Human qPCR Primer Pair (NM_004299)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22
Forward Sequence:	GCCATTGTAGGAGGTAGTGGGT
Reverse Sequence:	CGAAGGCTTTCCAGGCTCACAT
ACCN:	NM_004299 , NM_004299.1 , NM_004299.2 , NM_004299.3 , NM_004299.4 , BC006323 , BC006323.2 , BG722908 , BP228998 , BQ630948 , BT009918 , BX537833 , NM_004299.6
UniProt ID:	O75027
Synonyms:	ABC7; ASAT; Atm1p; EST140535
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)