

Product datasheet for **HP207660**

BNIP2 Human qPCR Primer Pair (NM_004330)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	663
Forward Sequence:	GGCTTAGACACACCGTCAGAGA
Reverse Sequence:	AACATACGCCAGCGTCGTCAT
ACCN:	<u>NM_004330</u> , <u>NM_004330.1</u> , <u>NM_004330.2</u> , <u>NM_004330.3</u> , <u>BC002461</u> , <u>BG257875</u> , <u>BX378770</u> , <u>NM_004330.4</u>
UniProt ID:	<u>Q12982</u>
Synonyms:	BNIP-2; NIP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)