

Product datasheet for **HP207661**

BNIP3L Human qPCR Primer Pair (NM_004331)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	665
Forward Sequence:	TGTGGAAATGCACACCAGCAGG
Reverse Sequence:	CTACTGGACCAGTCTGATACCC
ACCN:	<u>NM_004331</u> , <u>NM_004331.1</u> , <u>NM_004331.2</u> , <u>BC009603</u> , <u>BC009603.1</u> , <u>BC001559</u> , <u>BF337139</u> , <u>BG942896</u> , <u>BQ719112</u> , <u>BT019500</u> , <u>BT019501</u> , <u>NM_004331.3</u>
UniProt ID:	<u>O60238</u>
Synonyms:	BNIP3a; NIX
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)