

Product datasheet for **HP207663**

B Raf (BRAF) Human qPCR Primer Pair (NM_004333)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 673 |
| Forward Sequence: | AACGAGACCGATCCTCATCAGC |
| Reverse Sequence: | GGTAGCAGACAAACCTGTGGTTG |
| ACCN: | NM_004333 , NM_004333.1 , NM_004333.2 , NM_004333.3 , NM_004333.4 , BC101757 , BC101757.1 , BC038966 , BC112079 , BM263585 , NM_004333.6 |
| UniProt ID: | P15056 |
| Synonyms: | B-raf; B-RAF1; BRAF1; NS7; RAFB1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)