

Product datasheet for **HP207674**

Caspase 3 (CASP3) Human qPCR Primer Pair (NM_004346)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	836
Forward Sequence:	GGAAGCGAATCAATGGACTCTGG
Reverse Sequence:	GCATCGACATCTGTACCAGACC
ACCN:	NM_004346 , NM_004346.1 , NM_004346.2 , NM_004346.3 , BC016926 , BC016926.2 , BU753483 , BX647609 , NM_004346.4
UniProt ID:	P42574
Synonyms:	CPP32; CPP32B; SCA-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)