

Product datasheet for **HP207702**

CRF1 (CRHR1) Human qPCR Primer Pair (NM_004382)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1394
Forward Sequence:	TGCCAGGAGATCCTCAATGAGG
Reverse Sequence:	AGCCGCAGAAAGAGGACAAAGG
ACCN:	<u>NM_004382, NM_004382.1, NM_004382.3, NM_004382.4, BC096836, BC037967, NM_004382.5</u>
UniProt ID:	<u>P34998</u>
Synonyms:	CRF-R; CRF-R-1; CRF-R1; CRF1; CRFR-1; CRFR1; CRH-R-1; CRH-R1; CRHR; CRHR1L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)