

Product datasheet for **HP207704**

Versican (VCAN) Human qPCR Primer Pair (NM_004385)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1462
Forward Sequence:	TTGGACCTCAGGCGCTTTCTAC
Reverse Sequence:	GGATGACCAATTACTCAATCAC
ACCN:	<u>NM_004385</u> , <u>NM_004385.1</u> , <u>NM_004385.2</u> , <u>NM_004385.3</u> , <u>NM_004385.4</u> , <u>BC050524</u> , <u>BC172276</u> , <u>BP332680</u> , <u>BX641036</u> , <u>NM_004385.5</u>
UniProt ID:	<u>P13611</u>
Synonyms:	CSPG2; ERVR; GHAP; PG-M; WGN; WGN1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)