

## Product datasheet for **HP207709**

### Cathepsin H (CTSH) Human qPCR Primer Pair (NM\_004390)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1512
Forward Sequence:	TACCTTCGAGGTACTGGTCCCT
Reverse Sequence:	GGTGGAGAAAGTCCAGCAACTG
ACCN:	<a href="#">NM_004390</a> , <a href="#">NM_004390.1</a> , <a href="#">NM_004390.2</a> , <a href="#">NM_004390.3</a> , <a href="#">NM_004390.4</a> , <a href="#">BC002479</a> , <a href="#">BC002479.2</a> , <a href="#">BC018821</a> , <a href="#">BG775606</a> , <a href="#">BM849425</a> , <a href="#">NM_004390.5</a>
UniProt ID:	<a href="#">P09668</a>
Synonyms:	ACC-4; ACC-5; ACC4; ACC5; CPSB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)