

## Product datasheet for **HP207716**

### DDX6 Human qPCR Primer Pair (NM\_004397)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1656
Forward Sequence:	CCGAAATGGCTTATGCCGCAATC
Reverse Sequence:	GGAGATAGGTCTCTGCCAGCTT
ACCN:	<a href="#">NM_004397</a> , <a href="#">NM_004397.1</a> , <a href="#">NM_004397.2</a> , <a href="#">NM_004397.3</a> , <a href="#">NM_004397.4</a> , <a href="#">NM_004397.5</a> , <a href="#">BC065007</a> , <a href="#">BC026228</a> , <a href="#">BC039826</a> , <a href="#">BF810127</a> , <a href="#">NM_004397.6</a>
UniProt ID:	<a href="#">P26196</a>
Synonyms:	HLR2; IDDILF; P54; RCK; Rck/p54
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)