

Product datasheet for **HP207724**

Dnmt2 (TRDMT1) Human qPCR Primer Pair (NM_004412)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1787
Forward Sequence:	CACAAGTGGTGGCTGCCATTGA
Reverse Sequence:	GTCAAACCTCTCGAGTGAATGCC
ACCN:	<u>NM_004412</u> , <u>NM_004412.1</u> , <u>NM_004412.2</u> , <u>NM_004412.3</u> , <u>NM_004412.4</u> , <u>NM_004412.6</u> , <u>BC047733</u> , <u>BC047733.1</u> , <u>BU620523</u> , <u>BX537961</u> , <u>NM_004412.7</u>
UniProt ID:	<u>O14717</u>
Synonyms:	DMNT2; DNMT2; MHSAIIP; PUMET; RNMT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)