

Product datasheet for **HP207726**

Calcipressin 1 (RCAN1) Human qPCR Primer Pair (NM_004414)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1827
Forward Sequence:	GGAAACAAGTGGAAGATGCGACC
Reverse Sequence:	TGTCAGTCGCTGCGTGAATTC
ACCN:	NM_004414 , NM_004414.1 , NM_004414.2 , NM_004414.3 , NM_004414.4 , NM_004414.5 , NM_004414.6 , BC002864 , BC002864.2 , BT007363 , BX419260 , NM_004414.7
UniProt ID:	P53805
Synonyms:	ADAPT78; CSP1; DSC1; DSCR1; MCIP1; RCN1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)