

Product datasheet for **HP207737**

PHC1 Human qPCR Primer Pair (NM_004426)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1911
Forward Sequence:	TTCCTGCCAAGGCATCTCCAGT
Reverse Sequence:	AGTGTAGGAGCGGTACTGAAG
ACCN:	NM_004426 , NM_004426.1 , NM_004426.2 , BC002871 , BC017748 , BC073964 , BC157850 , BM473731 , BX640639 , NM_004426.3
UniProt ID:	P78364
Synonyms:	EDR1; HPH1; MCPH11; RAE28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)