

Product datasheet for **HP207741**

EGR3 Human qPCR Primer Pair (NM_004430)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1960
Forward Sequence:	GACTCGGTAGTCCATTACAATCAG
Reverse Sequence:	AGTAGGTCACGGTCTTGTTGCC
ACCN:	<u>NM_004430, NM_004430.1, NM_004430.2, BC104765, BC104765.1, BC112279, NM_004430.3</u>
UniProt ID:	<u>Q06889</u>
Synonyms:	EGR-3; PILOT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)