

## Product datasheet for **HP207749**

### Eph receptor A7 (EPHA7) Human qPCR Primer Pair (NM\_004440)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2045
Forward Sequence:	CTCGCAAGTGAGCGGAGTAATG
Reverse Sequence:	GTTGAGTAGGTCCGTTCCCTTTG
ACCN:	<a href="#">NM_004440</a> , <a href="#">NM_004440.1</a> , <a href="#">NM_004440.2</a> , <a href="#">NM_004440.3</a> , <a href="#">BC126125</a> , <a href="#">BC027940</a> , <a href="#">BC126151</a> , <a href="#">BC143857</a> , <a href="#">BC143858</a> , <a href="#">BM677016</a> , <a href="#">NM_004440.4</a>
UniProt ID:	<a href="#">Q15375</a>
Synonyms:	EHK-3; EHK3; EK11; HEK11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)