

Product datasheet for **HP207759**

Estrogen Related Receptor alpha (ESRRA) Human qPCR Primer Pair (NM_004451)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2101
Forward Sequence:	CCACTATGGTGTGGCATCCTGT
Reverse Sequence:	GGTGATCTCACACTCGTTGGAG
ACCN:	<u>NM_004451</u> , <u>NM_004451.1</u> , <u>NM_004451.2</u> , <u>NM_004451.3</u> , <u>NM_004451.4</u> , <u>BC063795</u> , <u>BC007915</u> , <u>BC011528</u> , <u>BC033701</u> , <u>BC073956</u> , <u>BC092470</u> , <u>BC093720</u> , <u>BC093722</u> , <u>NM_004451.5</u>
UniProt ID:	<u>P11474</u>
Synonyms:	ERR1; ERRa; ERRalpha; ESRL1; NR3B1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)