

Product datasheet for **HP207764**

EZH2 Human qPCR Primer Pair (NM_004456)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2146
Forward Sequence:	GACCTCTGTCTTACTTGTGGAGC
Reverse Sequence:	CGTCAGATGGTGCCAGCAATAG
ACCN:	<u>NM_004456</u> , <u>NM_004456.1</u> , <u>NM_004456.2</u> , <u>NM_004456.3</u> , <u>NM_004456.4</u> , <u>BC010858</u> , <u>BC010858.2</u> , <u>BG829016</u> , <u>BM974565</u> , <u>BU953788</u> , <u>NM_004456.5</u>
UniProt ID:	<u>Q15910</u>
Synonyms:	ENX-1; ENX1; EZH2b; KMT6; KMT6A; WVS; WVS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)