

Product datasheet for **HP207788**

GALNT3 Human qPCR Primer Pair (NM_004482)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2591
Forward Sequence:	GTTGCTAGGAGCAACAGTCGCA
Reverse Sequence:	CTGGACTTACGACAGCCGTGTA
ACCN:	<u>NM_004482</u> , <u>NM_004482.1</u> , <u>NM_004482.2</u> , <u>NM_004482.3</u> , <u>BC056246</u> , <u>BC113565</u> , <u>BC113567</u> , <u>BU567493</u> , <u>BX094158</u> , <u>BX647473</u> , <u>NM_004482.4</u>
UniProt ID:	<u>Q14435</u>
Synonyms:	GalNAc-T3; HFTC; HFTC1; HHS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)