

Product datasheet for **HP207794**

GPS2 Human qPCR Primer Pair (NM_004489)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2874
Forward Sequence:	GAGAAGCACCAGCTTTTCCTGC
Reverse Sequence:	GAACAGTCAGGCTCTGCTGGTA
ACCN:	NM_004489 , NM_004489.1 , NM_004489.2 , NM_004489.3 , NM_004489.4 , BC107738 , BC107738.1 , BC013652 , BC014192 , BC103901 , BC103902 , BC103903 , BT006998 , NM_004489.5
UniProt ID:	Q13227
Synonyms:	AMF-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)