

Product datasheet for **HP207797**

TFIIA2 (GTF2A2) Human qPCR Primer Pair (NM_004492)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2958
Forward Sequence:	AAATGCAGCACTGGCTCAGAGG
Reverse Sequence:	AAGTTCTGTCACCTCTCTGAATTC
ACCN:	<u>NM_004492</u> , <u>NM_004492.1</u> , <u>NM_004492.2</u> , <u>BC001919</u> , <u>BC001919.1</u> , <u>BC000287</u> , <u>BT007362</u> , <u>BX099306</u> , <u>NM_004492.3</u>
UniProt ID:	<u>P52657</u>
Synonyms:	HsT18745; T18745; TF2A2; TFIIA; TFIIA-12; TFIIA-gamma; TFIIAS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)