

## Product datasheet for **HP207806**

### HOXC6 Human qPCR Primer Pair (NM\_004503)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3223
Forward Sequence:	TTACCCCTGGATGCAGCGAATG
Reverse Sequence:	CCGCGTTAGGTAGCGATTGAAG
ACCN:	<a href="#">NM_004503</a> , <a href="#">NM_004503.1</a> , <a href="#">NM_004503.2</a> , <a href="#">NM_004503.3</a> , <a href="#">BC074845</a> , <a href="#">BC074845.2</a> , <a href="#">BC074844</a> , <a href="#">BF034547</a> , <a href="#">BF058345</a> , <a href="#">BF108967</a> , <a href="#">BF434131</a> , <a href="#">NM_004503.4</a>
UniProt ID:	<a href="#">P09630</a>
Synonyms:	CP25; HHO.C8; HOX3; HOX3C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)