

Product datasheet for **HP207808**

HSF2 Human qPCR Primer Pair (NM_004506)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3298
Forward Sequence:	GTTCAGCCAGTTTCGGAAGAGG
Reverse Sequence:	GAAGAAGCAGGGTCCCATCGA
ACCN:	NM_004506 , NM_004506.1 , NM_004506.2 , NM_004506.3 , BC112323 , BC112323.1 , BC121050 , BC005329 , BC064622 , BC121051 , BC128420 , NM_004506.4
UniProt ID:	Q03933
Synonyms:	HSF 2; HSTF 2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)