

Product datasheet for **HP207833**

NAP1L3 Human qPCR Primer Pair (NM_004538)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4675
Forward Sequence:	GAAATCCAGCCAAGAGTGGTTCC
Reverse Sequence:	TCCAGGATAGCATCTTCTCGTGG
ACCN:	NM_004538 , NM_004538.1 , NM_004538.2 , NM_004538.3 , NM_004538.4 , NM_004538.5 , BC034954 , BC034954.2 , BC094729 , NM_004538.6
UniProt ID:	Q99457
Synonyms:	MB20; NPL3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)